

GI Mentor Simulator

COVID-19 Disinfection Protocol

Note: When using the simulator, it is recommended to use protection equipment (disposable gloves and face mask).

Prior to Cleaning Simulator:

1. Power down the simulator and disconnect AC power.
2. Put on disposable gloves.


Cleaning Simulator Parts:

1. Dampen a new clean, microfiber cloth with a disinfectant (see table below). Ensure excess liquid is squeezed from the cloth. You may also use pre-dampened Isopropyl Alcohol wipes.
2. Gently wipe the moistened cloth/wipe on the surfaces to be cleaned.

Warning:

- Never clean a simulator while it is powered on or plugged in.
- Never spray/spill liquids of any kind directly onto the simulator.

Simulator Part	Disinfectant
1 Touchscreen Monitor	Standard IPA (Isopropyl Alcohol) at 70% alcohol concentration for cleaning and disinfecting. DO NOT use IPA at concentration >70%
2 Keyboard 3 Scopes 4 Simulator body	Standard IPA (Isopropyl Alcohol) for cleaning and disinfecting (at least 70%)



After Cleaning:

1. Make sure there is no visible moisture on the surfaces of the simulator before it is powered on or plugged in.
2. Remove and dispose gloves.
3. Wash your hands immediately after gloves are disposed.

Disclaimer: Though this disinfection protocol follows the guidance published by the [CDC](https://www.cdc.gov), 3D Systems cannot guarantee with 100% certainty that treatment per the above instructions will result in the total eradication of coronavirus.

Contact us for more information:

healthcare@3dsystems.com | www.simbionix.com

 **3D SYSTEMS**