

# PELVIC Mentor Simulator COVID-19 Disinfection Protocol

**When using the simulator, it is recommended to use protection equipment (disposable gloves and face mask).**

## Prior to Cleaning Simulator:

1. Power down the simulator and disconnect AC power.
2. Put on disposable gloves.

## Cleaning Simulator Parts:


1. Dampen a new clean microfiber cloth with the recommended disinfectant (see specifications below). Ensure excess liquid is squeezed from the cloth. You may also use pre-dampened isopropyl alcohol wipes.
2. Gently wipe the moistened surfaces to be cleaned.

### Warning:

- Never clean a simulator while it is powered on or plugged in.
- Never spray/spill liquids of any kind directly onto the simulator.

## PELVIC Mentor Simulator Parts and Disinfectants

Simulator Part	Disinfectant
<b>1</b> Computer Touchscreen	Standard IPA (Isopropyl Alcohol) at 70% alcohol concentration for cleaning and disinfecting. DO NOT use IPA at concentration >70%
<b>2</b> Pelvic Manikin - Front Vaginal Plate <b>3</b> Interchangeable Anatomical Parts (remove RF ID Tag prior to cleaning)	Warm water with soap for overall cleaning or Standard IPA for spot cleaning
<b>4</b> Pelvic Manikin - Plastic Frame <b>5</b> Platform Base <b>6</b> Control Unit	Standard IPA (Isopropyl Alcohol) for cleaning and disinfecting (at least 70%)




Contact us for more information:

healthcare@3dsystems.com | www.simbionix.com

## Additional Parts (refer to PELVIC Mentor Service Guide for disassembly)

Simulator Part	Disinfectant
<ol style="list-style-type: none"> <li>1 Finger Sensors</li> <li>2 Thimbles</li> <li>3 Speculum</li> <li>4 Pap Test Tools</li> <li>5 Switchbox</li> <li>6 Pelvic Skeleton</li> </ol>	<p>Standard IPA (Isopropyl Alcohol) for cleaning and disinfecting (at least 70%)</p>



## PELVIC Mentor Display Cart

Simulator Part	Disinfectant
<ol style="list-style-type: none"> <li>1 Computer Tray</li> <li>2 Keyboard and Mouse Tray</li> <li>3 Rotating Manikin Table</li> <li>4 Table Height Adjustment Mechanism</li> <li>5 Control Unit</li> </ol>	<p>Standard IPA (Isopropyl Alcohol) for cleaning and disinfecting (at least 70%)</p>



## After Cleaning:

1. Make sure there is no visible moisture on the surfaces of the simulator before it is powered on or plugged in.
2. Remove and dispose gloves.
3. Wash your hands immediately after gloves are disposed.

**Disclaimer:** Though this disinfection protocol follows the guidance published by the [CDC](https://www.cdc.gov), 3D Systems cannot guarantee with 100% certainty that treatment per the above instructions will result in the total eradication of coronavirus.

Contact us for more information:

healthcare@3dsystems.com | www.symbionix.com